

## PD-L1 antibody [MAB0871]

## Cat. No. GTX52447

Host	Rat
Clonality	Monoclonal
Isotype	lgG2
Applications	FCM, Neutralizing/Inhibition
Reactivity	Mouse

Package 200 μg

## Applications

## **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Assay dependent
Neutralizing/Inhibition	Assay dependent
Note: Pleaking the hinding of DD 11 and DD 1	

Note: Blocking the binding of PD-L1 and PD-1

Not tested in other applications.

Properties	
Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse PD-L1 transfected cells
Purification	Protein A/G purified From tissue culture supernatant
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 07 Page 1 of 1