

PDGF beta antibody [5F66]

Cat. No. GTX52467

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P
Reactivity	Human

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
IHC-P	Assay dependent

Not tested in other applications.

Calculated MW 27 kDa. ([Note](#))

Product Note This antibody shows approximately some degree (40%) cross-reactivity with PDGF-AB and no cross-reactivity with recombinant PDGF-AA.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	200 µg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human PDGF B recombinant protein
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX52467 WB Image

WB analysis of human placenta tissue lysates using GTX52467 PDGF beta antibody [5F66].


GTX52467 IHC-P Image

IHC-P analysis of formalin fixed human placenta tissue using GTX52467 PDGF beta antibody [5F66].



For full product information, images and publications, please visit our [website](#).