

CD62E antibody [9E7]

Cat. No. GTX52484

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2
Applications	WB, IHC-P, IP
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
IHC-P	1:25 - 1:50
IP	Assay dependent

Not tested in other applications.

Calculated MW	67 kDa. (Note)
Product Note	No cross reactivity to ICAM, and VCAM.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human CD62E.
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated


For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX52484 IHC-P Image

IHC-P analysis of formalin fixed human placenta tissue using GTX52484 CD62E antibody [9E7].



For full product information, images and publications, please visit our [website](#).