

Activin Receptor Type IIA antibody [6J16]

Cat. No. GTX52509

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2
Applications	WB, IHC-P
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:100 - 1:1000
IHC-P	1:50 - 1:100

Not tested in other applications.

Calculated MW 58 kDa. ([Note](#))

Product Note

This antibody was selected for its ability to detect human Activin Receptor type IIA. No cross reactivity was found to Activin RI and RIIb

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	200 µg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant Activin RIIA extra cellular domain
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated

For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 11 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX52509 IHC-P Image

IHC-P analysis of formalin fixed human placenta tissue using GTX52509 Activin Receptor Type IIA antibody [6J16].



For full product information, images and publications, please visit our [website](#).