

# NKp46 antibody [8F24]

**Cat. No. GTX52945**

|                     |                        |
|---------------------|------------------------|
| <b>Host</b>         | Mouse                  |
| <b>Clonality</b>    | Monoclonal             |
| <b>Isotype</b>      | IgG2                   |
| <b>Applications</b> | WB, IHC-P, IHC-Fr, FCM |
| <b>Reactivity</b>   | Human                  |

**Package**  
50 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500 - 1:1000       |
| IHC-P              | 1:25                 |
| IHC-Fr             | 1:20 - 1:100         |
| FCM                | 1:50 - 1:100         |

Not tested in other applications.

**Calculated MW** 34 kDa. ( [Note](#) )

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS                                                                                                                                                                                                                        |
| <b>Preservative</b>  | No preservatives                                                                                                                                                                                                           |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 500 µg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                 |
| <b>Immunogen</b>     | Human recombinant NKp46                                                                                                                                                                                                    |
| <b>Purification</b>  | Protein G purified<br>From tissue culture supernatant                                                                                                                                                                      |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |




For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

**DATA IMAGES**

**GTX52945 IHC-Fr Image**

IHC-Fr analysis of diseased human liver tissue using GTX52945 NKp46 antibody [8F24].


**GTX52945 IHC-P Image**

IHC-P analysis of human primary cutaneous melanoma using GTX52945 NKp46 antibody [8F24].



For full product information, images and publications, please visit our [website](#).