

PDGF Receptor alpha antibody [3G28]

Cat. No. GTX53096

Host	Rat
Clonality	Monoclonal
Isotype	IgG2
Applications	WB, IHC-P
Reactivity	Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:400 - 1:800
IHC-P	1:50 - 1:100

Not tested in other applications.

Calculated MW	123 kDa. (Note)
Product Note	This antibody recognizes mouse PDGF Receptor alpha in western blot. No cross-reactivity to mouse PDGF Receptor beta.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	200 µg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse PDGFR α extra-cellular domain
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX53096 IHC-P Image

IHC-P analysis of pancreas tissue sections from LPS exposed mouse using GTX53096 PDGF Receptor alpha antibody [3G28].



For full product information, images and publications, please visit our [website](#).