

VEGFB antibody [6E72]

Cat. No. GTX53099

Host	Rat
Clonality	Monoclonal
Isotype	lgG2
Applications	WB, IHC-P
Reactivity	Mouse

Package 100 μg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:400 - 1:1000
IHC-P	1:50 - 1:100
Not tested in other applications.	

 Calculated MW
 22 kDa. (Note)

 Product Note
 This antibody recognizes mouse VEGF-B 186 but not 167.

Properties	
Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	$200 \ \mu g/ml$ (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse VEGF-B 186 recombinant protein
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 07 Page 1 of 2



DATA IMAGES



GTX53099 IHC-P Image

 $IHC-P\ analysis\ of\ lung\ tissue\ sections\ from\ LPS\ exposed\ mouse\ using\ GTX53099\ VEGFB\ antibody\ [6E72].$

Left : Primary antibody Right : Negative control



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 07 Page 2 of 2