

VEGFB antibody [6E72]

Cat. No. GTX53099

Host	Rat
Clonality	Monoclonal
Isotype	IgG2
Application	WB, IHC-P
Reactivity	Mouse

Package
100 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:400 - 1:1000
IHC-P	1:50 - 1:100

Not tested in other applications.

Calculated MW	22 kDa. (Note)
Product Note	This antibody recognizes mouse VEGF-B 186 but not 167.

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	200 µg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse VEGF-B 186 recombinant protein
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX53099 IHC-P Image

IHC-P analysis of lung tissue sections from LPS exposed mouse using GTX53099 VEGFB antibody [6E72].

Left : Primary antibody

Right : Negative control



For full product information, images and publications, please visit our [website](#).