

HAI-2 antibody [RM0089-7F52]

Cat. No. GTX53179

Host	Rat
Clonality	Monoclonal
Isotype	lgG2
Applications	WB
Reactivity	Mouse

Package 100 μg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:1000
Not tested in other applications.	

Calculated MW 28 kDa. (<u>Note</u>)

Properties	
Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	$500 \mu g/ml$ (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse recombinant of the common extracellular region of HAI-2 isoforms
Purification	Protein A/G purified From tissue culture supernatant
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our website.

Date 2025 / 11 / 03 Page 1 of 1