

TIM-3 antibody [6F46]

Cat. No. GTX53225

Host	Rat
Clonality	Monoclonal
Isotype	IgG2
Applications	WB
Reactivity	Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:1000

Not tested in other applications.

Calculated MW 31 kDa. ([Note](#))

Product Note This antibody detects mouse of TIM-3 in western blotting. This antibody does not cross-react with mouse TIM-1 and TIM-2.

Properties

Form Liquid

Buffer PBS

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 500 µg/ml (Please refer to the vial label for the specific concentration.)

Immunogen Mouse recombinant protein of TIM-3

Purification Protein A/G purified
From tissue culture supernatant

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).