

## Epiregulin antibody [6H12]

Cat. No. GTX53324

|              |            |
|--------------|------------|
| Host         | Rat        |
| Clonality    | Monoclonal |
| Isotype      | IgG2       |
| Applications | ELISA      |
| Reactivity   | Mouse      |

Package  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| ELISA              | 1:1000 - 1:2000      |

Not tested in other applications.

## Properties

|               |                                                                                                                                                                                                                                                                                                                                             |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                                                                                                                                      |
| Buffer        | PBS                                                                                                                                                                                                                                                                                                                                         |
| Preservative  | No preservatives                                                                                                                                                                                                                                                                                                                            |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.                                                                                                                  |
| Concentration | 200 µg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                                                                                                                                  |
| Immunogen     | Mouse recombinant protein Epiregulin                                                                                                                                                                                                                                                                                                        |
| Purification  | Protein G purified<br>From tissue culture supernatant                                                                                                                                                                                                                                                                                       |
| Conjugation   | Unconjugated                                                                                                                                                                                                                                                                                                                                |
| Note          | <p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p> |



For full product information, images and publications, please visit our [website](#).