

ID1 antibody

Cat. No. GTX53624

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human, Mouse, Rat

Package
400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
ICC/IF	1:10-1:50
IHC-P	1:10-1:50
FCM	1:25

Note : Permeabilization : 90% methanol for 10 min

Not tested in other applications.

Calculated MW 16 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 66-93 amino acids from the Central region of human ID1.
Purification	Protein A purified, followed by peptide affinity purification.
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX53624 IHC-P Image

IHC-P analysis of human pancreas tissue using GTX53624 ID1 antibody.



GTX53624 WB Image

WB analysis of U251 cell lysate (35ug/lane) using GTX53624 ID1 antibody.



GTX53624 WB Image

WB analysis of mouse heart tissue lysate (35ug/lane) using GTX53624 ID1 antibody.

Dilution : 1:1000



For full product information, images and publications, please visit our [website](#).