

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody Titer Assay Kit (Mouse)

Cat. No. GTX536582

Applications ELISA

Package
96 test

PRODUCT

Summary

The GeneTex SARS-CoV-2 (COVID-19) Nucleocapsid Antibody Titer Assay Kit (Mouse) is an in vitro tool that measures the mouse antibody titer against SARS-CoV-2 nucleocapsid protein with the lowest detectable dose concentration of standard nucleocapsid antibody being 6.8 pg/ml (mouse monoclonal antibody).

Applications

Application Note

The suitable sample types include but not limit to purified mouse antibodies and serum. The suggested sample loading amount is 50µl per well. Please see the protocol for detailed procedure.

Properties


Preservative 0.05% ProClin 300

Storage The entire kit can be stored at 2-8°C for up to 6 months.

Conjugation Horseradish peroxidase(HRP)

Note For *In vitro* laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES



GTX536582 ELISA Image

ELISA detection of mouse monoclonal SARS-CoV-2 (COVID-19) nucleocapsid antibody (GTX635808) using SARS-CoV-2 (COVID-19) Nucleocapsid Antibody Titer Assay Kit (Mouse) (GTX536582).



For full product information, images and publications, please visit our [website](#).