

## CPI-17 antibody [4H10]

Cat. No. GTX53770

|              |                |
|--------------|----------------|
| Host         | Mouse          |
| Clonality    | Monoclonal     |
| Isotype      | IgG2a          |
| Applications | WB, FCM, ELISA |
| Reactivity   | Human          |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500                |
| FCM                | Assay dependent      |
| ELISA              | Assay dependent      |

Not tested in other applications.

Calculated MW 17 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS                                                                                                                                                                                                                        |
| Preservative  | 0.1% Sodium azide                                                                                                                                                                                                          |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| Immunogen     | Recombinant human CPI-17 (1-147aa) purified from E. coli                                                                                                                                                                   |
| Purification  | By protein-G affinity chromatography                                                                                                                                                                                       |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX53770 WB Image

WB analysis of A549 lysate (30ug) using CPI-17 antibody [4H10] at a dilution of 1:500.



For full product information, images and publications, please visit our [website](#).