

# SARS-CoV-2 Neutralizing Antibody ELISA Kit (Omicron BA.4 / BA.5 / BF.7 / BQ.1 / XBB.1.5)

**Cat. No. GTX538288**

**Applications** Neutralizing /Inhibition

**Reactivity** SARS Coronavirus 2

References ( 2 )

Package

96 test

## PRODUCT

### Summary

The GeneTex SARS-CoV-2 Neutralizing Antibody ELISA Kit (Omicron BA.4/5 / BF.7 / BQ.1 / XBB1.5) is an in vitro assay for qualitative SARS-CoV-2 NAb screening. The kit includes recombinant His-tagged SARS-CoV-2 Spike RBD proteins for the Omicron variant BA.4/5, Omicron variant BF.7, Omicron variant BQ.1, and Omicron variant XBB 1.5, for Omicron variant-specific neutralization analysis.

### Properties


**Storage** The entire kit can be stored at 2-8°C for up to 6 months.

**Conjugation** Horseradish peroxidase(HRP)

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


**Note** Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



### GTX538288 Neutralizing /Inhibition Image

Inhibition analysis of SARS-CoV-2 (COVID-19) Spike RBD antibody [HL1004] (GTX635793) (104-13333 pM) using SARS-CoV-2 Neutralizing Antibody ELISA Kit (Omicron BA.4 / BA.5 / BF.7 / BQ.1 / XBB.1.5).



### GTX538288 Neutralizing /Inhibition Image

Inhibition analysis of SARS-CoV-2 (COVID-19) Spike RBD Omicron antibody [HL1866] (GTX637591) (104-13333 pM) using SARS-CoV-2 Neutralizing Antibody ELISA Kit (Omicron BA.4 / BA.5 / BF.7 / BQ.1 / XBB.1.5).



For full product information, images and publications, please visit our [website](#).



**GTX538288 Neutralizing /Inhibition Image**



For full product information, images and publications, please visit our [website](#).