

SARS-CoV-2 Neutralizing Antibody ELISA Kit (Omicron BA.4 / BA.5 / BF.7 / BQ.1 / XBB.1.5)

Cat. No. GTX538288

Applications Neutralizing/Inhibition

Reactivity SARS Coronavirus 2

References (2)

Package

96 test

PRODUCT

Summary

The GeneTex SARS-CoV-2 Neutralizing Antibody ELISA Kit (Omicron BA.4/5 / BF.7 / BQ.1 / XBB1.5) is an in vitro assay for qualitative SARS-CoV-2 NAb screening. The kit includes recombinant His-tagged SARS-CoV-2 Spike RBD proteins for the Omicron variant BA.4/5, Omicron variant BF.7, Omicron variant BQ.1, and Omicron variant XBB 1.5, for Omicron variant-specific neutralization analysis.

Properties

Storage The entire kit can be stored at 2-8°C for up to 6 months.


Conjugation Horseradish peroxidase(HRP)

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX538288 Neutralizing/Inhibition Image

Inhibition analysis of SARS-CoV-2 (COVID-19) Spike RBD antibody [HL1004] (GTX635793) (104-13333 pM) using SARS-CoV-2 Neutralizing Antibody ELISA Kit (Omicron BA.4 / BA.5 / BF.7 / BQ.1 / XBB.1.5).



GTX538288 Neutralizing/Inhibition Image

Inhibition analysis of SARS-CoV-2 (COVID-19) Spike RBD Omicron antibody [HL1866] (GTX637591) (104-13333 pM) using SARS-CoV-2 Neutralizing Antibody ELISA Kit (Omicron BA.4 / BA.5 / BF.7 / BQ.1 / XBB.1.5).



For full product information, images and publications, please visit our [website](https://www.genetex.com).



GTX538288 Neutralizing/Inhibition Image



For full product information, images and publications, please visit our [website](#).