

IL33 (pro) antibody

Cat. No. GTX54067

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 31 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	IL-33 antibody was raised against a 19 amino acid synthetic peptide from near the amino terminus of human IL-33. The immunogen is located within amino acids 50 - 100 of IL-33.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

DATA IMAGES



GTX54067 WB Image

WB analysis of human lymph node tissue lysate using GTX54067 IL33 (pro) antibody.
Working concentration : (A) 1 and (B) 2 μ g/ml



GTX54067 IHC-P Image

IHC-P analysis of human lymph node tissue using GTX54067 IL33 (pro) antibody.
Working concentration : 20 μ g/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).

Date 2026 / 10 / 06 Page 2 of 2