

GDF15 antibody

Cat. No. GTX54079

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IP
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
IHC-P	1:50 - 1:100
IP	1:50 - 1:200

Not tested in other applications.

Calculated MW 34 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 49-308 of human GDF15 (NP_004855.2).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX54079 IP Image

IP analysis of HT-1080 cell lysate using GTX54079 GDF15 antibody.

Antibody amount : 3 μ g / 200 μ g lysate


Dilution : 1:1000


GTX54079 WB Image

WB analysis of various sample lysates using GTX54079 GDF15 antibody.

Dilution : 1:1000

Loading : 25 μ g per lane


GTX54079 IHC-P Image

IHC-P analysis of human placenta tissue using GTX54079 GDF15 antibody.

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).