

Respiratory Syncytial virus Nucleoprotein ELISA pair [HL1248 / HL1246]

Cat. No. GTX541617

Package
1 pair

PRODUCT

Product Content

Cat No.	Type	Host	Conjugation	Package
GTX636648	Detection Antibody	Rabbit	Unconjugated	100µl
GTX636650	Capture Antibody	Rabbit	Unconjugated	100µl

Properties

Form Liquid

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration batch dependent (Please refer to the vial label for the specific concentration.)

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 15 Page 1 of 1