

Nucleostemin antibody

Cat. No. GTX54163

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IP
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	1:10 - 1:100
IP	1:50 - 1:100

Not tested in other applications.

Calculated MW 62 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human GNL3 (NP_996562.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX54163 ICC/IF Image

ICC/IF analysis of MCF7 cells using GTX54163 Nucleostemin antibody.
Blue : DAPI


GTX54163 WB Image

WB analysis of 293T cell lysate using GTX54163 Nucleostemin antibody.
Dilution : 1:1000
Loading : 25µg per lane


GTX54163 IP Image

IP analysis of 293T cell lysate using GTX54163 Nucleostemin antibody.
Antibody amount : 3µg / 200µg lysate
Dilution : 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).