

CD36 antibody [CRF D2712]

Cat. No. GTX54400

Host	Mouse
Clonality	Monoclonal
Isotype	IgA
Applications	WB, ICC/IF, FCM, IP, Functional Assay
Reactivity	Mouse, Rat

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
FCM	Assay dependent
IP	Assay dependent
Functional Assay	Assay dependent

Not tested in other applications.

Calculated MW 53 kDa. ([Note](#))

Product Note Monoclonal antibody CRF D-2712 blocks CD36 binding. It is not known whether the antibody is activating or neutralizing.

Properties

Form	Liquid
Buffer	Tissue culture supernatant
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Purification	Unpurified
Endotoxin	Low endotoxin
Conjugation	Unconjugated

For full product information, images and publications, please visit our [website](#).

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).