

L-FABP antibody [K5A6]

Cat. No. GTX54433

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-Fr, ELISA
Reactivity	Human, Baboon

Package
50 µg

PRODUCT

Summary The monoclonal antibody K5A6 is useful to detect ischemic areas of human liver.

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

WB	Assay dependent
----	-----------------

IHC-Fr	Assay dependent
--------	-----------------

ELISA	Assay dependent
-------	-----------------

Not tested in other applications.

Calculated MW 14 kDa. ([Note](#))

Properties

Form Liquid

Buffer Filter-sterilized PBS, 0.1% BSA

Preservative 0.02% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.

Concentration 100 µg/ml (Please refer to the vial label for the specific concentration.)

Purification Purified IgG1

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).