

C5 / C5a antibody [561]

Cat. No. GTX54438

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ELISA, Functional Assay
Reactivity	Human

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Functional Assay	Assay dependent

Note : The antibody 561 can be used as detector in ELISA

Not tested in other applications.

Product Note Monoclonal antibody 561 recognizes the N-terminus of C5. The antibody reacts both with intact C5 (190 kDa) as with C5a (115 kDa).

Properties

Form	Liquid
Buffer	Filter-sterilized PBS, 0.1% BSA
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	100 µg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human C5/C5a
Purification	Purified IgG2a
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).