

Myeloperoxidase antibody [2D4]

Cat. No. GTX54497

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	FACS, ELISA
Reactivity	Mouse, Rat

Package
50 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FACS	Assay dependent
ELISA	Assay dependent

**Note : Cells are fixed with 4% paraformaldehyde and permeabilized with 0.1% saponine.
The antibody can be used as a detection antibody.**

Not tested in other applications.

PROPERTIES

Form	Liquid
Buffer	Filter-sterilized PBS, 0.1% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	100 µg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse myeloperoxidase
Purification	Protein G purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).