

## CD55 antibody [RDIII-7]

Cat. No. GTX54502

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	WB, ICC/IF, IHC-Fr, FACS, Neutralizing/Inhibition
Reactivity	Rat

Package  
50 µg

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-Fr	Assay dependent
FACS	Assay dependent
Neutralizing/Inhibition	Assay dependent

**Note : RDIII-7 antibody completely blocked the protective effect of expressed rat DAF.**

Not tested in other applications.

## PROPERTIES

Form	Liquid
Buffer	Filter-sterilized PBS, 0.1% BSA
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	100 µg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	NIH-3T3 cells expressing GPI-anchored rat DAF
Purification	Purified IgG1
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).