

## CD62E antibody [ENA1]

Cat. No. GTX54505

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	IHC-Fr, IP, ELISA, Functional Assay
Reactivity	Human, Baboon, Chimpanzee, Cynomolgus monkey, Rhesus Monkey

Package  
50 µg

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-Fr	Assay dependent
IP	Assay dependent
ELISA	Assay dependent
Functional Assay	Assay dependent

Not tested in other applications.

**Calculated MW** 67 kDa. ( [Note](#) )**Product Note** The antibody inhibits the adhesion of granulocytes both neutrophilic and eosinophilic.

## PROPERTIES

Form	Liquid
Buffer	Filter-sterilized PBS, 0.1% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	100 µg/ml (Please refer to the vial label for the specific concentration.)
Purification	Purified IgG1
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).