

NFkB p100 (phospho Ser866) antibody

Cat. No. GTX54550

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	1:50 - 1:100
IHC-P	1:50 - 1:100

Not tested in other applications.

Calculated MW 97 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A phospho specific peptide corresponding to residues surrounding S866 of human NFkB2
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX54550 WB Image

WB analysis of MDA-MB435 cell treated by TNF alpha using GTX54550 NFκB p100 (phospho Ser866) antibody.

Loading : 25μg per lane


GTX54550 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTX54550 NFκB p100 (phospho Ser866) antibody.


GTX54550 ICC/IF Image

ICC/IF analysis of methanol-fixed HeLa cells using GTX54550 NFκB p100 (phospho Ser866) antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).