

PDGF beta antibody

Cat. No. GTX54575

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

References (2)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	Assay dependent
IHC-P	1:50 - 1:200

Not tested in other applications.

Calculated MW 27 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.05% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-241 of human PDGFB (NP_002599.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX54575 WB Image

WB analysis of various cell lysates using GTX54575 PDGF beta antibody.

Dilution : 1:1000


Loading : 25µg per lane


GTX54575 IHC-P Image

IHC-P analysis of human placenta tissue using GTX54575 PDGF beta antibody.

Dilution : 1:150

Antigen retrieval : 10 mM citrate buffer, pH 6.0


GTX54575 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX54575 PDGF beta antibody.

Dilution : 1:150

Blue : DAPI



For full product information, images and publications, please visit our [website](https://www.genetex.com).