

Staphylococcus aureus PBP2a (Methicillin-resistant) overexpression Expi293F whole cell lysate

Cat. No. GTX548834

Species Staphylococcus aureus

Package

200 µg

PRODUCT

Summary

Whole-cell lysate of Expi293F cells overexpressing Staphylococcus aureus PBP2a (Methicillin-resistant) protein (full-length, C-terminal 3x Flag-tagged). The translated amino acid sequence is based on ALJ10987.1. Lysate prepared in PBS with 1% Tween-20 buffer with no reducing agents and is unboiled. Please see matched control Expi293F lysate ([GTX548835](#))

Properties

Form Liquid

Buffer PBS, 1% Tween-20

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration batch dependent (Please refer to the vial label for the specific concentration.)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).