

Adenosine A2b Receptor antibody

Cat. No. GTX54903

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:1000
ICC/IF	1:100 - 1:500
IHC-P	1:100 - 1:200

Not tested in other applications.

Calculated MW 36 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.42% Potassium Phosphate, 0.87% NaCl, 30% Glycerol
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of Adenosine A2b Receptor. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX54903 ICC/IF Image

ICC/IF analysis of formalin-fixed SH-SY5Y cells using GTX54903 Adenosine A2b Receptor antibody.

Red : Primary antibody

Blue : DAPI


Permeabilization : 0.1% Triton X-100 in TBS for 5-10 minutes



GTX54903 IHC-P Image

IHC-P analysis of formalin fixed human pancreas tissue section using GTX54903 Adenosine A2b Receptor antibody.

Antigen retrieval : Heat mediated antigen retrieval with sodium citrate buffer (pH 6.0)



GTX54903 WB Image

WB analysis of SH-SY5Y (A), NIH3T3 (B), rat kidney (C) whole cell lysates using GTX54903 Adenosine A2b Receptor antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).