

SHC1 (phospho Tyr239) antibody

Cat. No. GTX54986

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:1000
IHC-P	1:100 - 1:200

Not tested in other applications.

Calculated MW 63 kDa. ([Note](#))

Product Note

This antibody is raised against human SHC1 phosphorylated at Tyr349 of p66Shc (Uniprot #P29353-1). Based on sequence homology, it is corresponding to human SHC1 phosphorylated at Tyr239 of p52Shc (Uniprot #P29353-2).

Properties

Form	Liquid
Buffer	0.42% Potassium Phosphate, 0.87% NaCl, 30% Glycerol
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of SHC1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 08 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX54986 WB Image

WB analysis of EGF-treated HeLa (A), EGF-treated HUVEC(B) whole cell lysates using GTX54986 SHC1 (phospho Tyr239) antibody.



GTX54986 IHC-P Image

IHC-P analysis of formalin fixed human prostate cancer tissue section using GTX54986 SHC1 (phospho Tyr239) antibody.

Antigen retrieval : Heat mediated antigen retrieval with sodium citrate buffer (pH 6.0)



For full product information, images and publications, please visit our [website](https://www.genetex.com).