

c-Abl (phospho Tyr204) antibody

Cat. No. GTX54998

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:1000
IHC-P	1:100 - 1:200

Not tested in other applications.

Calculated MW 123 kDa. ([Note](#))**Product Note** This antibody is raised against human c-Abl (isoform b) phosphorylated at Tyr204. Based on sequence homology, it is predicted to react with c-Abl (isoform a) when phosphorylated at the corresponding residues - c-Abl (isoform a) Tyr185.

Properties

Form	Liquid
Buffer	0.42% Potassium Phosphate, 0.87% NaCl, 30% Glycerol
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of c-Abl. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTx54998 WB Image

WB analysis of colchicine-treated HeLa (A), mouse kidney (B), colchicine-treated PC12 (C) whole cell lysates using GTx54998 c-Abl (phospho Tyr204) antibody.



GTx54998 IHC-P Image

IHC-P analysis of formalin fixed human breast cancer tissue section using GTx54998 c-Abl (phospho Tyr204) antibody.

Antigen retrieval : Heat mediated antigen retrieval with sodium citrate buffer (pH 6.0)



For full product information, images and publications, please visit our [website](https://www.genetex.com).