

Chk1 (phospho Ser296) antibody

Cat. No. GTX55053

| | |
|---------------------|------------|
| Host | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Applications | WB, IHC-P |
| Reactivity | Human |

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB | 1:500 - 1:1000 |
| IHC-P | 1:100 - 1:200 |

Not tested in other applications.

Calculated MW 54 kDa. ([Note](#))


Properties

| | |
|----------------------|--|
| Form | Liquid |
| Buffer | 0.42% Potassium Phosphate, 0.87% NaCl, 30% Glycerol |
| Preservative | 0.01% Sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.) |
| Immunogen | KLH-conjugated synthetic peptide encompassing a sequence within the center region of Chk1. The exact sequence is proprietary. |
| Purification | Purified by antigen-affinity chromatography |
| Conjugation | Unconjugated |
| Note | For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product. |




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

**GTX55053 IHC-P Image**

IHC-P analysis of formalin fixed human breast cancer tissue section using GTX55053 Chk1 (phospho Ser296) antibody.

Antigen retrieval : Heat mediated antigen retrieval with sodium citrate buffer (pH 6.0)

**GTX55053 WB Image**

WB analysis of UV-treated HeLa (A) whole cell lysates using GTX55053 Chk1 (phospho Ser296) antibody.



For full product information, images and publications, please visit our [website](#).