

## TCF3 / E2A (phospho Thr355) antibody

Cat. No. GTX55358

|              |            |
|--------------|------------|
| Host         | Rabbit     |
| Clonality    | Polyclonal |
| Isotype      | IgG        |
| Applications | IHC-P      |
| Reactivity   | Human      |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| IHC-P              | 1:50-1:100           |

Not tested in other applications.


## Properties

|               |                                                                                                                                                                                                                                                                                                                                             |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                                                                                                                                      |
| Buffer        | PBS, 150mM NaCl, 50% Glycerol                                                                                                                                                                                                                                                                                                               |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                                                                                                                                          |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.                                                                                                                  |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                                                                                                                                    |
| Immunogen     | Peptide sequence around phosphorylation site of threonine 355 (P-S-T(p)-P-V) derived from human TCF3.                                                                                                                                                                                                                                       |
| Purification  | Purified by antigen-affinity chromatography. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.                                                                                                                                                                                                       |
| Conjugation   | Unconjugated                                                                                                                                                                                                                                                                                                                                |
| Note          | <p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p> |



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX55358 IHC-P Image

IHC-P analysis of human brain tissue using GTX55358 TCF3 / E2A (phospho Thr355) antibody.

Left : Primary antibody

Right : Primary antibody pre-incubated with the antigen specific peptide



For full product information, images and publications, please visit our [website](#).