

Heme Oxygenase 1 antibody

Cat. No. GTX55656

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	1:50 - 1:200
IHC-P	1:50 - 1:200

Not tested in other applications.

Calculated MW 33 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein of human HMOX1
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX55656 WB Image

WB analysis of various sample lysates using GTX55656 Heme Oxygenase 1 antibody.

Dilution : 1:1000


Loading : 25µg per lane


GTX55656 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX55656 Heme Oxygenase 1 antibody.

Blue : DAPI

Dilution : 1:100


GTX55656 IHC-P Image

IHC-P analysis of human liver tissue using GTX55656 Heme Oxygenase 1 antibody.

Dilution : 1:200



For full product information, images and publications, please visit our [website](https://www.genetex.com).