

## TGF beta Receptor 2 antibody

Cat. No. GTX55814

|              |                   |
|--------------|-------------------|
| Host         | Rabbit            |
| Clonality    | Polyclonal        |
| Isotype      | IgG               |
| Applications | WB, ICC/IF, IHC-P |
| Reactivity   | Human, Mouse, Rat |

References ( 3 )

Package

100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500 - 1:2000       |
| ICC/IF             | 1:50 - 1:200         |
| IHC-P              | 1:50 - 1:200         |

Not tested in other applications.

Calculated MW 65 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 50% Glycerol                                                                                                                                                                                                          |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | Recombinant fusion protein containing a sequence corresponding to amino acids 23-280 of human TGFBR2 (NP_003233.4).                                                                                                        |
| Purification  | Purified by affinity chromatography                                                                                                                                                                                        |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTx55814 IHC-P Image

IHC-P analysis of human appendix tissue using GTx55814 TGF beta Receptor 2 antibody.  
Dilution : 1:100




### GTx55814 WB Image

WB analysis of various sample lysates using GTx55814 TGF beta Receptor 2 antibody.  
Dilution : 1:1000  
Loading : 25µg per lane



### GTx55814 IHC-P Image

IHC-P analysis of human breast cancer tissue using GTx55814 TGF beta Receptor 2 antibody.  
Dilution : 1:100



### GTx55814 ICC/IF Image

ICC/IF analysis of NIH/3T3 cells using GTx55814 TGF beta Receptor 2 antibody.  
Blue : DAPI  
Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).