

## Visfatin antibody

Cat. No. GTX55840

|              |                   |
|--------------|-------------------|
| Host         | Rabbit            |
| Clonality    | Polyclonal        |
| Isotype      | IgG               |
| Applications | WB, ICC/IF, IHC-P |
| Reactivity   | Human, Mouse, Rat |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500 - 1:2000       |
| ICC/IF             | 1:50 - 1:200         |
| IHC-P              | 1:50 - 1:200         |

Not tested in other applications.

Calculated MW 56 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 50% Glycerol                                                                                                                                                                                                          |
| Preservative  | 0.05% ProClin 300                                                                                                                                                                                                          |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human NAMPT (NP_005737.1).                                                                                                          |
| Purification  | Purified by affinity chromatography                                                                                                                                                                                        |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTx55840 WB Image

WB analysis of various sample lysates using GTx55840 Visfatin antibody.


Dilution : 1:1000

Loading : 25µg per lane



### GTx55840 ICC/IF Image


ICC/IF analysis of HeLa cells using GTx55840 Visfatin antibody.



### GTx55840 IHC-P Image

IHC-P analysis of mouse spleen tissue using GTx55840 Visfatin antibody.

Dilution : 1:100



### GTx55840 IHC-P Image

IHC-P analysis of rat ovary tissue using GTx55840 Visfatin antibody.

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).