

ZEB1 antibody

Cat. No. GTX55847

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

References (2)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:1000
IHC-P	1:50 - 1:200

Not tested in other applications.

Calculated MW 124 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein of human ZEB1
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX55847 WB Image

WB analysis of Jurkat cell lysate using GTX55847 ZEB1 antibody.

Dilution : 1:1000


Loading : 25µg per lane



GTX55847 IHC-P Image

IHC-P analysis of human gastric cancer tissue using GTX55847 ZEB1 antibody.

Dilution : 1:200




GTX55847 WB Image

WB analysis of normal (control) and knockout (KO) HeLa cell lysate using GTX55847 ZEB1 antibody.

Dilution : 1:1000

Loading : 25µg per lane



GTX55847 WB Image

The data was published in the journal J Exp Clin Cancer Res in 2019. [PMID: 31661013](https://pubmed.ncbi.nlm.nih.gov/31661013/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).