

## PI3 kinase C2 gamma antibody

Cat. No. GTX55876

|              |                   |
|--------------|-------------------|
| Host         | Rabbit            |
| Clonality    | Polyclonal        |
| Isotype      | IgG               |
| Applications | WB                |
| Reactivity   | Human, Mouse, Rat |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500 - 1:1000       |

Not tested in other applications.

Calculated MW 166 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | 0.42% Potassium Phosphate, 0.87% NaCl, 30% Glycerol                                                                                                                                                                        |
| Preservative  | 0.01% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of PI3 kinase C2 gamma. The exact sequence is proprietary.                                                                               |
| Purification  | Purified by antigen-affinity chromatography                                                                                                                                                                                |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX55876 WB Image

WB analysis of CEM (A), HeLa (B), NIH3T3 (C), PC12 (D) whole cell lysates using GTX55876 PI3 kinase C2 gamma antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).