

# APR2 antibody

**Cat. No. GTX56097**

|                     |                   |
|---------------------|-------------------|
| <b>Host</b>         | Rabbit            |
| <b>Clonality</b>    | Polyclonal        |
| <b>Isotype</b>      | IgG               |
| <b>Applications</b> | WB                |
| <b>Reactivity</b>   | Human, Mouse, Rat |

**Package**  
100 µl

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500 - 1:1000       |

Not tested in other applications.


## Properties

|                      |                                                                                                                                                                                                                                                                                                                                      |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                                                                                                                               |
| <b>Buffer</b>        | 0.42% Potassium Phosphate, 0.87% NaCl, 30% Glycerol                                                                                                                                                                                                                                                                                  |
| <b>Preservative</b>  | 0.01% Sodium azide                                                                                                                                                                                                                                                                                                                   |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.                                                                                                           |
| <b>Concentration</b> | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                                                                                                                                     |
| <b>Immunogen</b>     | KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of APR2. The exact sequence is proprietary.                                                                                                                                                                                                        |
| <b>Purification</b>  | Purified by antigen-affinity chromatography                                                                                                                                                                                                                                                                                          |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                                                                                                                                         |
| <b>Note</b>          | For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.<br><br>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product. |



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX56097 WB Image

WB analysis of HEK293T (A), NS-1 (B), H9C2 (C) whole cell lysates using GTX56097 APR2 antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).