

JNK (phospho Tyr185) antibody

Cat. No. GTX57151

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µl

Applications

Application Note

WB: 1:500-1:1000

Product Note

This antibody is raised against JNK1/2 phosphorylated at Tyr185. Based on sequence homology, it is predicted to react with JNK3 when phosphorylated at the corresponding residues - Tyr223.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic phosphopeptide derived from human JNK1 around the phosphorylation site of Tyrosine 185.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).