

TGF beta Receptor I antibody

Cat. No. GTX57156

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µl

Applications

Application Note

WB: 1:500 - 1:1000, ICC/IF: 1:50 - 1:200

Calculated MW 56 kDa. ([Note](#))


Product Note TGFβ RI (E161) pAb detects endogenous levels of TGFβ RI protein.

Properties


Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide, corresponding to amino acids 150-200 of Human TGFβ RI.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx57156 WB Image

Western blot (WB) analysis of TGFβ RI pAb at 1:500 dilution Lane1:Hela cell lysate Lane2:LO2 cell lysate
Lane3:Raw264.7 cell lysate Lane4:Mouse kidney tissue lysate Lane5:Rat heart tissue lysate


GTx57156 WB Image

Western blot (WB) analysis of TGFβ RI pAb at 1:1000 dilution Lane1:The Liver tissue lysate of Rat(40ug)
Lane2:The Liver tissue lysate of mouse(40ug) Lane3:LO2 whole cell lysate(40ug) Lane4:The Heart tissue lysate
of Rat(40ug)



For full product information, images and
publications, please visit our [website](https://www.genetex.com).