

CD57 antibody [IHC539]

Cat. No. GTX57181

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

Package
500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100-1:200

Note : Heat-induced epitope retrieval (HIER) at pH 9 for 10 - 30 minutes

Not tested in other applications.

Properties

Form	Liquid
Buffer	Tris, 1% BSA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant Human CD57 protein.
Purification	Protein A purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX57181 IHC-P Image

IHC-P analysis of human tonsil tissue using GTX57181 CD57 antibody [IHC539]



For full product information, images and publications, please visit our [website](#).