

IDH1 (R132H Mutant Specific) antibody [IHC132]

Cat. No. GTX57185

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-P
Reactivity	Human

References (4)

Package

500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100-1:200

Note : Heat-induced epitope retrieval (HIER) at pH 9 for 10 - 30 minutes

Not tested in other applications.

Properties

Form	Liquid
Buffer	Tris, 1% BSA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide against the C-Terminus of IDH1, with a point mutation at residual 132 where an Arginine R is changed to Histidine. protein.
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX57185 IHC-P Image

IHC-P analysis of human astrocytoma tissue using GTX57185 IDH1 (R132H Mutant Specific) antibody [IHC132]



For full product information, images and publications, please visit our [website](#).