

TFE3 antibody [IHC672]

Cat. No. GTX57200

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-P
Reactivity	Human

Package
500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100-1:200

Note : Heat-induced epitope retrieval (HIER) at pH 9 for 10 - 30 minutes

Not tested in other applications.

Properties

Form	Liquid
Buffer	Tris, 1% BSA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant Human TFE3 protein.
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](https://www.genetex.com).

DATA IMAGES



GTX57200 IHC-P Image

IHC-P analysis of human kidney tissue using GTX57200 TFE3 antibody [IHC672]



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 2 of 2