

Human HSFY1 protein, His tag

Cat. No. GTX57258-pro

Species Human**Package**

100 µg

Applications

Observed MW (kDa) 45 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 0.4M Urea, 10% Glycerol**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMahVSSE TQDVSPKDEL TASEASTRSP
LCEHTFPGDS DLRSMIIEHA FQVLSQGSLL ESPSYTVCVS EPDKDDDFLS LNFPRKLWKI VESDQFKSIS WDENGTCIVI NEELFKKEIL
ETKAPYRIFQ TDAIKSFVRQ LNLYGFSKIQ QNFQRS AFLA TFLSEEKES VLSKLKYYN PNFKRGYPQL LVRVKRRIGV KNASPISTLF
NEDFNKKHFR AGANMENHNS ALAAEASEES LFSASKNLNM PLTRESSVRQ IIANSSVPIR SGFPPSPST SVGPSEQIAT
DQHAILNQLT TIHMHSHTY MQARGHIVNF ITTTTSQYHI ISPLQNGYFG LTVEPSAVPT RYPLVSVNEA PYRNMLPAGN
PWLQMPTIAD RSAAPHSRLA LQPSPLDKYH PNYN

Expression System E. coli**Purity** > 80% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX57258-pro Image

3μg Human HSFY1 protein (GTX57258-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).