

Human Fank1 protein, His tag

Cat. No. GTX57275-pro

Species Human

Package
100 µg

Applications

Observed MW (kDa) 38 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 150mM NaCl, 20% Glycerol, 1mM DTT, 0.1mM PMSF

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMEPQKIM PPSKPHPPV GKVTHHSIEL YWDLEKKAKR QGPQEQWFRF SIEEDPKMH TYGIITGYA TKHVVEGLEP RTLYRFRKLV TSPSGECEYS PLVSVSTTRE PISSEHLHRA VSVNDEDLV RILQGGRVKV DVPNKFGFTA LMVAAQKGYT RLVKILVSNG TDVNLKNGSG KDSLMLACYA GHLDVVKYLR RHGASWQARD LGGCTALHWA ADGGHCSVIE WMIKDGCEVD VDTGSGWTP LMRVSAVSGN QRVASLLIDA GANVNVKDRN GKTPLMVAVL NNHEELVQLL LDKGADASVK NEFGKGVLEM ARVFDRQSVV SLEERKKKQ RPKKSCVC

Expression System E. coli

Purity > 85% by SDS-PAGE.


Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx57275-pro Image

3µg Human Fank1 protein (GTx57275-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).