

Human OSR2 protein, His tag

Cat. No. GTX57302-pro
Species Human

Package
100 µg

Applications
Observed MW (kDa) 36 kDa.

Properties
Form Liquid

Buffer 20mM Tris-HCl, 0.4M Urea, 10% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMGSKALP APIPLHPSLQ LTNYSFLQAV NTFPATVDHL QGLYGLSAVQ TMHMHNHWTLG YPNVHEITRS TITEMAAAQG LVDARFPFPA LPFTTHLFHP KQGAIAHVLP ALHKDRPRFD FANLAVAATQ EDPPKMGDLS KLSPGLGSPi SGLSKLTPDR KPSRGRLPSK TKKEFICKFC GRHFTKSYNL LIHERHTHDE RPYTCDICHK AFRRQDHLRD HRYIHSKEKP FKQCQCGKGF CQSRTLAVHK TLHMQTSSPT AASSAAKCSG ETVICGGTA

Expression System E. coli


Purity > 80% by SDS-PAGE.

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX57302-pro Image

3μg Human OSR2 protein (GTX57302-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).