

Human TOM1L2 protein, His tag

Cat. No. GTX57358-pro

Species	Human
---------	-------

Package

50 µg

Applications

Observed MW (kDa)	56 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	PBS, 10% Glycerol, 1mM DTT
--------	----------------------------

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSEFMEFLL GNPFPSTPVGQ CLEKATDGSL QSEDWTLNME ICDIINETEE GPKDAIRALK KRLNGNRNYR EVMLALTVE TCVKNCGRHF HILVANRDFI DSVLVKIISP KNNPPTIVQD KVLALIQAWA DAFRSPDLT GVWHIYEELK RKGVEFPMAD LDALSPIHTP QRSVPEVDPA ATMPRSQSQQ RTSAGSYSSP PPAPYSAPQA PALSVTGPIT ANSEQIARLR SELDVVRGNT KVMSEMLTEM VPGQEDSSDL ELLQELNRTC RAMQQRIVEL ISRVSNEEVT EELLHVNDL NNVFLRYERF ERYRSGRSVQ NASNGVLNEV TEDNLIDLGP GSPAVVSPMV GNTAPPSSLS SQLAGLDLGT ESVSGTLSSL QQCNPDRGFD MFAQTRGNSL AEQRKTVTYE DPQAVGGLAS ALDNRKQSSE GIPVAQPSVM DDIEVWLRTD LKGGDLEEGV TSEEFDKFLE ERAKAAEMVP DLPSPPEAP APASNPSGRK KPERSEDALF AL

Expression System	E. coli
-------------------	---------

Purity	> 80% by SDS-PAGE.
--------	--------------------

Conjugation	Unconjugated
-------------	--------------

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX57358-pro Image

3μg Human TOM1L2 protein (GTX57358-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).