

## Human Spred1 protein, His tag

Cat. No. GTX57371-pro

|         |       |
|---------|-------|
| Species | Human |
|---------|-------|

|         |       |
|---------|-------|
| Package | 50 µg |
|---------|-------|

## Applications

|                   |         |
|-------------------|---------|
| Observed MW (kDa) | 50 kDa. |
|-------------------|---------|

## Properties

|      |        |
|------|--------|
| Form | Liquid |
|------|--------|

|        |                             |
|--------|-----------------------------|
| Buffer | 20mM Tris-HCl, 10% Glycerol |
|--------|-----------------------------|

|              |                  |
|--------------|------------------|
| Preservative | No preservatives |
|--------------|------------------|

|         |                                                                                                                                                                                                                            |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|               |                                                                            |
|---------------|----------------------------------------------------------------------------|
| Concentration | 0.5 mg/ml (Please refer to the vial label for the specific concentration.) |
|---------------|----------------------------------------------------------------------------|

## Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMSSEETAT SDNDNSYARV RAVVMTRDDS SGGWLPLGGS GLSSVTVFKV PHQEENGCAD FFIGERLRD KVVVLECMK KDLYNKVTP TFHHWKIDDK KFGTLFQSPA DARAFDRGIR RAIEDISQGC PESKNEAEGA DDLQANEEDS SSSLVKDHLE QQETVVTSEP YRSSNIRPSP FEDLNARRVY MQSQANQITF GQPGLDIQR SMEYVQRQIS KECGSLKSN RVPLKSIRHV SFQDEDEIVR INPRDILIRR YADYRHPDMW KNDLERDDAD SSIQFSKPSD KKS DYLYSCG DETKLSSPKD SVVFKTQPSS LKIKKSKRRK EDGERSRCVY CQERFNHEEN VRGKCQDAPD PIKRCIQVS CMLCAESMLY HCMSDSEGDF SDPCSCDTS DDKFCLRWLAL VALSFIVPCM CCYVPLRMCH RCGEACGCCG GKHKAAAG

|                   |         |
|-------------------|---------|
| Expression System | E. coli |
|-------------------|---------|

|        |                    |
|--------|--------------------|
| Purity | > 80% by SDS-PAGE. |
|--------|--------------------|

|             |              |
|-------------|--------------|
| Conjugation | Unconjugated |
|-------------|--------------|

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX57371-pro Image

3μg Human Spred1 protein (GTX57371-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).